POST OP SHOE

SUP3125

To see all the FAQs,
visit coretechortho.link/u4a

coretech®

QUICK START GUIDE

INTENDED USE

Provides post-surgical support and protection for foot and toe recovery following bunionectomy, hammer toe correction, metatarsal fractures, or diabetic wound management. Rigid sole stabilizes foot and shifts weight away from surgical site to promote healing.

HOW TO WEAR

1. **Slide foot into shoe** – ensure heel is positioned correctly
2. **Secure bottom strap** across top of foot
3. **Fasten ankle strap** around ankle
4. **Adjust for comfort** – straps should be snug but not restrict circulation

NOTE: Shoe works on left or right foot. Can be worn with or without socks.

STRAP ADJUSTMENT

Straps are designed to fit over bulky dressings. If straps are excessively long, you may trim them with scissors. Cut small amounts at a time and round the corners to prevent fraying. Do not cut near the attachment points.

CARE INSTRUCTIONS

- Spot clean with damp cloth and mild soap.
- Do NOT machine wash or tumble dry.
- Allow to air dry completely before wearing.
- Insole can be removed and replaced if needed.

WARNINGS ⚠

- CONSULT YOUR PHYSICIAN before use if you have circulation problems
- NOT FOR HIGH-IMPACT ACTIVITIES – walking only
- SINGLE PATIENT USE ONLY – do not share
- DISCONTINUE USE if pain increases or redness/swelling occurs
- DO NOT OVERTIGHTEN straps – may restrict blood flow

WARRANTY

1 year limited warranty

coretech®

by vivehealth

Distributed by

vive health
8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

060126